

ICT+

isolator
cleaning tool

MAXIMUM CARE FOR MINI ENVIRONMENTS



 **SATURIX**[®]
by FG CLEAN WIPES



ISOLATOR CLEANING TOOL+

A HIGHER LEVEL OF SAFETY AND EFFICIENCY

Conventional isolator cleaning tools require extensive handling to attach and remove pads, which takes time and poses a significant risk for cross contamination. Saturix ICT+, with patent-pending designs, not only eliminates the need to change pads by hand, but also provides a complete system for used-pad management with superior overall usability.

It will help any compounding pharmacy be more compliant with USP standards and ensure the clean and sterile preparation of drugs for the safety of patients and employees.

TOUCH-FREE PAD APPLICATION AND REMOVAL

- Clean pads click on the frame, dirty pads click off the frame—all touch-free
- Reduces cross contamination risk while increasing efficiency

ERGONOMIC, EASY-TO-USE HYBRID DESIGN

- Asymmetrical frame maximizes coverage and effectively cleans tight spaces
- Offset yoke on the frame extends reach—easy to clean high areas

99.9% PARTICULATE RESIDUE REMOVAL

- Microfiber cleaning pads remove 99.9% of particles down to .1 micron
- Mechanically removes disinfectant residues and other harmful particulates from the surface



Touch-free pad attachment and removal



Precision dosing of pads in package




Used-pad storage within the hood

HYBRID FRAME DESIGN
Easy to clean high and tight spaces

OPTIMIZED SURFACE COVERAGE WITH MINIMIZED CARBON FOOTPRINT

Each pad covers 138 ft² (12.8 m²)

 100% recyclable pad backer further reduces waste (30% of total replacement pad weight)

Removes **99.9%** of Particulate Residue Down to .1 micron*

* Based on 3rd party lab testing

SATURIX ICT+ ISOLATOR CLEANING TOOL

ICT+ Hardware

| SKU# | Product Description | Size | Qty/Pack | Packs/Case |
|--------|---|--|----------|------------|
| IFH19A | Medical Grade Stainless Steel Frame & Handle Kit, Auto-Clavable | Handle: 19" (480mm) Frame: 8.2" x 2.6" (208 mm x 66 mm) | 1 | 1 |
| IFH19 | Aluminum Frame & Stainless Steel Handle Kit | Handle: 19" (480mm) Frame: 8.2" x 2.6" (208 mm x 66 mm) | 1 | 1 |

ICT+ Microfiber Pads

| | | | | |
|--------|---|------------------------------|---|----|
| IPDP8S | Validated Sterile Microfiber Pad in Package | 8.2" x 2.6" (208 mm x 66 mm) | 4 | 15 |
| IPDB8S | Validated Sterile Microfiber Pads | 8.2" x 2.6" (208 mm x 66 mm) | 2 | 30 |
| IPDP8 | Non-Sterile Microfiber Pads in Package | 8.2" x 2.6" (208 mm x 66 mm) | 4 | 15 |
| IPDB8 | Non-Sterile Microfiber Pads | 8.2" x 2.6" (208 mm x 66 mm) | 4 | 15 |
| IPA30S | Validated Sterile 70% Isopropyl Alcohol | 30 ml, 1 FL OZ | 4 | 15 |



ICT+ isolator
cleaning tool

SATURIX[®]
by FG CLEAN WIPES

2255 Westover Road, Chicopee, MA 01022
+1 (413) 598-8344 • saturix.com

A Filtration Group Company

ICT-001.2